

CHROMOSPACE ● ● ● ● ● ● ● ● ● ●



● ● ● ● ● ● ● ● ● ●  
**MOSTLEDS**

PHOTOMEDECINE®, THE SCIENCE OF ACTIVE LIGHT

**CHROMOSPACE®**  
THE LARGEST SPACE OF  
PHOTOBIMODULATION

The CHROMOSPACE® integrates a unique meshing of more than 4 LED/cm<sup>2</sup>. OSRAM LED assure a homogenous and optimal use for 50 000 hours.

With the CHROMOSPACE®, LED phototherapy enters in a large space of efficiency.

The CHROMOSPACE® is the only LED light generator offering the most successful device of light therapy. The CHROMOSPACE® enables face and body treatments at the same time, to optimize and shorten the time of LED sessions on large treatment areas.

## TECHNICAL SPECIFICATIONS

NUMBER OF LED	22 400
TREATMENT AREA	8 532,8 cm <sup>2</sup>
DIMENSIONS	H 130 cm x L 220 cm x l 80 cm
5 ACTIVE COLORS	BLUE (460 nm) GREEN (525 nm) YELLOW (595 nm) RED (645 nm) NEAR INFRARED (850 nm)
TREATMENT PROTOCOLS	Pre-set protocols Well-being programs Expert mode: to create your own protocols
CERTIFICATION	CE, ISO 13 485 Quality Insurance for medical devices
MADE IN FRANCE -2 YEARS WARRANTY	

## BENEFITS

### PREMIUM IMAGE

Offer to your patients the latest generation of LED technology and an exceptional dedicated space to LED light.

### EXCEPTIONAL PROFITABILITY

Thanks to the high number of treated patients in a minimum of time.

### COMMUNICATION PACK

We provide communication tools to inform your patients (leaflets, posters, Internet SEO, videos...).

### LAUNCHING ASSISTANCE

We can provide a personalized support: tariff advices, open-house, VIP evening, commercial assistance, etc...

### QUALIFY TRAINING

- Initial training
- Advance training
- Continuous training at Paris pilot centre

### SIMPLE USE

- Hand-free device (no dedicated staff)
- No consumables (LED = 50 000 working hours)
- Energy savings

The CHROMOSPACE® device is exclusively reserved to medical and paramedical professions, clinics, medicalized and thermal spas.

[www.mostleds.com](http://www.mostleds.com)

2 rue de Lisbonne - 75008 Paris - Tel : +(33)1 42 93 67 82

